


# MVD-Enferplex™ Bovine BCV/BRSV

01D16 (465 Tests)

Test for the *in vitro* detection of antibodies to *Bovine Coronavirus* and *Bovine Respiratory Syncytial virus*.

**For *in vitro* veterinary diagnostic and research use only**



Manufactured & Distributed by:  
Enfer Scientific,  
Unit T, M7 Business Park,  
Newhall, Naas,  
Co. Kildare,  
Ireland.  
Tel.: +353 459 83800

---

## Contents Page

1.0	General Information .....	3
2.0	Intended Use.....	3
3.0	Principle of the Procedure .....	3
4.0	Reagents.....	3
5.0	Materials and Equipment required but not provided .....	4
6.0	Warnings and Precautions.....	5
7.0	Health and Safety Information.....	5
8.0	Preparation of Reagents.....	6
9.0	Sample & Control Preparation .....	7
10.0	Test Protocol.....	7
11.0	Results.....	8
12.0	Limitations of the Procedure .....	8
13.0	Disclaimer and Reservation of Rights.....	8
14.0	Recommended Plate Layout .....	9

## 1.0 General Information

### ***Bovine Coronavirus***

*Bovine coronaviruses* (BCoV) are members of the Coronaviridae family in the Nidovirus order, subgroup 2a, which are enveloped, single stranded, positive sense RNA viruses with shaped projections on the surface, (4). BCoV infects the upper and lower respiratory tract and intestine. It sheds in both feces and nasal secretions (5). Bovine CoV is involved in the etiology of 3 distinct clinical syndromes in cattle: calf diarrhea, winter dysentery with hemorrhagic diarrhea in adults, and respiratory infections in cattle of various ages including the bovine respiratory disease complex or shipping fever of feed lot cattle, (3). Economic losses are the primary concern of newborn or adult BCoV infection, as the morbidity rates are high and may reach up to 100% of the herd. In dairy herds suffering a BCoV outbreak milk production may not return to normal for several weeks, or even during that lactation period, resulting in significant losses for the milk industry (6).

### ***Bovine Respiratory Syncytial virus***

*Bovine Respiratory syncytial virus* (BRSV) is an enveloped, non-segmented, negative stranded RNA virus. This virus is classified within the genus Pneumovirus of the family Paramyxoviridae, (1). The virus is transmitted by aerosols and infects respiratory tract mucosal cells. The clinical signs of BRSV infection are associated with pulmonary localisation of the virus which is mainly in the cranio-ventral lobes of the lung, (2). BRSV infection is widespread and has major economic impact and is regarded as one of the most important causes of respiratory tract disease, especially in young calves. An infection can cause respiratory distress, fever, anorexia, pulmonary and subcutaneous emphysema and can lead to secondary bacterial pneumonia and death, (3).

## 2.0 Intended Use

The MVD-Enferplex™ Bovine BCV/BRSV assay is a semi-quantitative luminescent (emission) immunological method for the detection of various antibodies in bovine serum, individual and bulk milk. The MVD-Enferplex™ Bovine BCV/BRSV kit is intended for *in vitro* veterinary diagnostic use and research purposes only.

## 3.0 Principle of the Procedure

The MVD-Enferplex™ Bovine BCV/BRSV assay is a semi-quantitative enzyme immunoassay based on the sequential addition of bovine serum, individual or bulk milk to a multiple antigen coated plate, followed by antibody-enzyme conjugate and a chemiluminescent substrate.

Upon incubation of the test sample in the multiple antigen coated well, antibodies specific to the various diseases form complexes with the immobilized antigens.

This step is followed by a wash step with 1X Wash Buffer solution, and then Sheep anti-bovine sera labelled with HRPO (Horseradish peroxidase) is added, forming an antigen-antibody-conjugate-peroxidase complex. Next, unbound conjugate is washed away and a chemiluminescent substrate is used to generate the light signal and the image is captured. The image is analysed and data reduced to determine sample status in the appropriate Enferplex™ Macro.

## 4.0 Reagents

Reagent pack 01D16 contains sufficient material for 465 tests. The reagent pack is stored at 2-8°C. Note the storage conditions for individual components.

	Component	Function	Quantity	Storage Requirement
1.	BCV/BRSV Antibody Capture Plate	Plate used for antibody capture	5 x 96-well plates	2-8°C in sealed foil pouch
2.	20X Wash Buffer	Washes off unbound antibody/reagent	1 x 500ml of 20X concentrate solution (contains 0.05% Proclin 300® & 0.02% Protectol BN®)	2-8°C
3.	BMD Sample Diluent	Buffer for dilution of serum, or milk	2 x 500ml working strength solution (contains 0.05% Proclin 300® & 0.02% Protectol BN®)	2-8°C
4.	BMD Conjugate Diluent	Conjugate diluent for dilution of the BMD Concentrate Conjugate	1 x 500ml of working strength solution (contains 0.05% Proclin 300® & 0.02% Protectol BN®)	2-8°C
5.	Concentrate Conjugate	Undiluted Sheep anti-bovine IgG antibody	1 x microvial containing 0.1ml	2-8°C
6.	Multi-Lite A	Chemiluminescent substrate for peroxidase when combined with Multi-Lite B	1 x bottle containing 15ml solution	+15°C to 30°C (Room Temperature, RT)
7.	Multi-Lite B	Chemiluminescent substrate for peroxidase when combined with Multi-Lite A	1 x bottle containing 15ml solution	+15°C to 30°C
8.	BCV/BRSV Negative Control	Non-reactive with antibody capture plate, used as a control	1 x microvial containing 0.1ml	-15°C -to -25°C
9.	BCV/BRSV Positive Control	Non-infectious, Reactive with antibody capture plate, used as a control	1 x microvial containing 0.1ml	-15°C -to -25°C

## 5.0 Materials and Equipment required but not provided

- Microplate Incubator/shaker thermostated at 37°C ± 2°C and capable of shaking at 900rpm
- Quansys Biosciences Q-View Imager
- Device for the delivery and aspiration of wash solution
- High quality deionised, distilled, or reverse osmosis water
- Microplate Cover Seals
- Precision Micropipettes and Multichannel micropipettes of appropriate volume and disposable tips
- Reagent dispensing trays

- Glass containers for dilution of the concentrate conjugate
- Glass or polypropylene containers for dilution of other reagents
- Polypropylene tubes/plates for dilution of the samples

## 6.0 Warnings and Precautions

- 6.1 The reagents are solely for *in vitro* veterinary diagnostic and research only use on bovine serum, bulk and individual milk samples. For professional use only.
- 6.2 Please refer to the manufacturer's safety data sheets and the product labelling for information on potentially hazardous components.
- 6.3 Do not perform the test in the presence of reactive vapours (acids, alkalis, aldehydes) or dust, which could alter the enzymatic activity of the conjugate.
- 6.4 Use perfectly washed glassware, rinsed in distilled/deionised water or preferably disposable material.
- 6.5 Use a new pipette tip for each sample.
- 6.6 Do not modify the test procedure or substitute reagents from other manufacturers.
- 6.7 Do not use the reagents beyond the stated expiry date and do not intermix components from different kit lots. Microbiological contamination of reagents must be avoided as this may reduce the life of the product and cause erroneous results.
- 6.8 Use separate dispensing trays for each reagent used in the assay. Avoid cross contamination of the Multi-Lite solution with the diluted conjugate solution. Do not pour unused Multi-Lite solution back into the Multi-Lite bottles.
- 6.9 Do not allow plates to sit for more than 3 minutes between wash steps and the addition of reagents.
- 6.10 Do not expose the substrate to strong light or oxidizing agents.
- 6.11 Allow the reagents to adjust to room temperature (RT) (+18°C to 30°C), for approximately 30 minutes before use. Immediately after use, return all reagents to their appropriate storage conditions.
- 6.12 All reagents must be prepared in either clean glass, or polypropylene bottles. Care must be taken to avoid cross contamination of reagents.
- 6.13 Washing of the wells is an essential step of the procedure; ensure that all wells are completely filled, then completely emptied. Do not adjust the recommended washing steps. Inadequate washing can give incorrect results.

## 7.0 Health and Safety Information

The 20X Wash Buffer must be handled with care. Please note hazard identified on individual container label.

20X Wash Buffer contains 2-methyl-2H-isothiazol-3-one, which is classified as per EC Directive EC 1272/2008 Skin. Sens 1 – H317. The following are the appropriate Hazard (H) and Precautionary (P) Statements.



- H317 May cause an allergic skin reaction
- P261 Avoid breathing dust/fume/gas/mist/vapours/spray.
- P272 Contaminated work clothing should not be allowed out of the workplace
- P302 & P352 IF ON SKIN: Wash with plenty of soap and water
- P321 Specific treatment (see on this label)
- P333 & P313 If skin irritation or rash occurs: Get medical advice/attention.
- P363 Wash contaminated clothing before reuse.

Safety Data Sheets are available upon request.

## 8.0 Preparation of Reagents

Component	Method	Storage of Prepared Reagents
<b>Working Strength Wash Buffer (1X)</b> Requires: 1. 20X Wash Buffer 2. Deionised or distilled water	<ol style="list-style-type: none"> <li>Dilute 1 part 20X Wash Buffer in 19 parts deionised or distilled water.</li> <li>Mix thoroughly, e.g. For 4L of working strength 1X Wash Buffer, add 200ml 20X Wash Buffer to 3800ml of water</li> </ol>	1 month at RT or at 2-8°C
<b>Working Strength Conjugate</b> Requires: 1. BMD Conjugate Diluent 2. Concentrate Conjugate	<ol style="list-style-type: none"> <li>Prepare only the required volume for the number of tests to be carried out.</li> <li>400µl of working strength conjugate is required for 8 wells.</li> <li>5ml of working strength conjugate is required for 1 plate.</li> <li>Dilute the Concentrate Conjugate 1/25000 in BMD Conjugate Diluent, e.g. to 75ml of BMD Conjugate Diluent, add 3µl of Concentrate Conjugate.</li> <li>Mix by inversion. Invert a minimum of 8 times</li> </ol>	Prepare 15-20mins before use
<b>Multi-Lite Solution</b> Requires: 1. Multi-Lite A 2. Multi-Lite B	<ol style="list-style-type: none"> <li>Prepare only the required volume for the number of tests to be carried out.</li> <li>400µl of Multi-Lite Solution is required for 8 wells. 5ml of Multi-Lite Solution is required for 1 plate.</li> <li>Add equal parts of Multi-Lite A to Multi-Lite B in either a clean glass or plastic disposable vessel. For example, add 1ml of Multi-Lite A to 1ml of Multi-Lite B.</li> <li>Mix by inversion twice.</li> </ol>	Prepare 15-20mins before use and store in the dark
<b>Working Strength BCV/BRSV Negative Control</b> Requires: 1. BCV/BRSV Negative Control 2. BMD Sample Diluent	<ol style="list-style-type: none"> <li>Prepare only the required volume for the number of tests to be carried out.</li> <li>50µl of working strength BCV/BRSV Negative Control is required per plate.</li> <li>Dilute BCV/BRSV Negative Control 1/300 in the BMD Sample Diluent, e.g. to 1.2ml of BMD Sample Diluent, add 4µl of BCV/BRSV Negative Control.</li> <li>Mix by inversion a minimum of 2 times.</li> </ol>	Store at 2-8°C and use within 8hrs of preparation
<b>Working Strength BCV/BRSV Positive Control</b> Requires: 1. BCV/BRSV Positive Control 2. BMD Sample Diluent	<ol style="list-style-type: none"> <li>Prepare only the required volume for the number of tests to be carried out.</li> <li>50µl of working strength BCV/BRSV Positive Control is required per plate.</li> <li>Dilute BCV/BRSV Positive Control 1/300 in the BMD Sample Diluent, e.g. to 1.2ml of BMD Sample Diluent, add 4µl of BCV/BRSV Positive Control.</li> <li>Mix By inversion a minimum of 2 times.</li> </ol>	Store at 2-8°C and use within 8hrs of preparation

## 9.0 Sample & Control Preparation

Bring all specimens to room temperature prior to testing. All samples and controls must be added to the BCV/BRSV antibody capture plate at approximately the same time, therefore use of a transfer/master plate is recommended to add samples and controls to first, and then transfer to the BCV/BRSV antibody capture plate.

### 9.1 Controls

- 9.1.1 Add 50µl of BMD Sample Diluent to A1 of the test plate.
- 9.1.2 Add 50µl of BCV/BRSV Negative Control to B1 of the test plate.
- 9.1.3 Add 50µl of BCV/BRSV Positive Control to C1 of the test plate.

### 9.2 Serum or Plasma Samples

- 9.2.1 Fresh, refrigerated, or previously frozen serum can be tested. Icteric, lipemic, haemolysed, heat treated and contaminated sera may cause erroneous results.
- 9.2.2 If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage periods greater than 24 hours, freeze the serum at -20°C or below.
- 9.2.3 Specimens containing precipitate may yield inconsistent test results and such specimens must be clarified prior to testing.
- 9.2.4 Ensure the serum sample is mixed thoroughly before addition to the BMD Sample Diluent.
- 9.2.5 The samples are prepared to a 1:300 dilution by adding for example 4µl of the serum to 1.2ml of the BMD Sample Diluent.
- 9.2.6 Mix the prepared sample by inversion or by pipetting up and down a minimum of 2 times.

### 9.3 Bulk and individual Milk Samples

- 9.3.1 Whole milk samples can be used after centrifugation for 15 minutes at 2000 x g or left to stand if refrigerated (2-8°C). No pre-treatment is needed for defatted milk.
- 9.3.2 If specimens are not immediately tested, they should be refrigerated at 2-8°C. For storage periods greater than 24 hours, freeze the milk at -20°C or below.
- 9.3.3 Bulk milk samples are prepared to a 1/5 dilution by adding for example 18µl of the bulk milk to 90µl of the BMD Sample Diluent.
- 9.3.4 Mix the prepared sample by inversion or by pipetting up and down a minimum of 2 times.
- 9.3.5 Individual milk samples are prepared to a 1/10 dilution by adding for example 10µl of the individual milk to 100µl of the BMD Sample Diluent.
- 9.3.6 Mix the prepared sample by inversion or by pipetting up and down a minimum of 2 times.

## 10.0 Test Protocol

- 10.1 All samples are tested in singlicate. Refer to 13.0 for recommended plate layout.
- 10.2 Remove the BCV/BRSV antibody capture microplate from the protective packaging. For each single plate, transfer 50µl of the controls and samples into the wells of the BCV/BRSV antibody capture plate.
- 10.3 Cover the microplate with a microplate cover seal.
- 10.4 Incubate the microplate, shaking, for 60 minutes at 37 ± 2°C.
- 10.5 Remove the microplate cover seal, and wash the wells 6 times with 200/250µl of 1X wash buffer.
- 10.6 Dry by inversion on absorbent paper.
- 10.7 Add 50µl of working strength conjugate to each well. Cover with a microplate cover seal.
- 10.8 Incubate the microplate, shaking, for 60 minutes at 37 ± 2°C.

- 10.9** Remove the microplate cover seal, and wash the wells 6 times with 200/250µl of 1X wash buffer.
- 10.10** Dry by inversion on absorbent paper.
- 10.11** Add 50µl of the substrate solution to each well of the microplate. Immediately read the plate on the Q-View Imager set at 45 seconds exposure time.

## **11.0 Results**

### **11.1 Validation of Test Performance**

Each plate must be considered separately when calculating and interpreting results of the assay. The control results must be validated before the sample results can be interpreted. The criteria for the Buffer, BCV/BRSV Negative Control and BCV/BRSV Positive Control are all contained within the appropriate 'Enferplex™ BCV/BRSV Macro' provided and the results are calculated automatically.

### **11.2 Acceptable Range of Control Results**

If the criteria for the controls are not met, the assay is invalid and must be repeated.

### **11.3 Interpretation of Results**

#### **Negative Result**

Samples giving a 'Negative' result in the macro are considered non-reactive in the BCV/BRSV assay.

#### **Positive Result**

Samples giving a 'Positive' result in the macro are considered reactive in the BCV/BRSV assay.

## **12.0 Limitations of the Procedure**

Enfer Scientific complies with the quality system standard ISO9001:2008.

As with any biological test, this test may give a false positive or a false negative result owing to local conditions. A test should be interpreted in the context of all available clinical, historical, and epidemiological information relevant to the animal(s) under test.

A negative result with a semi-qualitative immunological method does not preclude the possibility of infection with *Bovine Coronavirus* or *Bovine Respiratory Syncytial Virus*.

Any change or modification of the procedure might affect the results.

*Responsibility for test interpretation and consequent animal husbandry decisions rests solely with the user and any consulting veterinarian and appropriate animal health advisors or authorities. Enfer Scientific accepts no responsibility for any loss or damage, howsoever caused, arising out of the interpretation of test results.*

## **13.0 Disclaimer and Reservation of Rights**

Enfer Scientific gives no warranty of any kind, whether expressed or implied, in regard to the carrying out of the MVD-Enferplex™ Bovine BCV/BRSV assay or for the stability and storage of the MVD-Enferplex™ BCV/BRSV kit, or for the procedure used. Without prejudice to the foregoing, Enfer Scientific disclaims all responsibility for merchantability and fitness for use after it leaves Enfer Scientific. Enfer Scientific shall not be liable, under any circumstances, for damages, direct or consequential.



## 14.0 Recommended Plate Layout

	1	2	3	4	5	6	7	8	9	10	11	12
<b>A</b>	<b>Buffer</b>	<b>S6</b>	<b>S14</b>	<b>S22</b>	<b>S30</b>	<b>S38</b>	<b>S46</b>	<b>S54</b>	<b>S62</b>	<b>S70</b>	<b>S78</b>	<b>S86</b>
<b>B</b>	<b>NC</b>	<b>S7</b>	<b>S15</b>	<b>S23</b>	<b>S31</b>	<b>S39</b>	<b>S47</b>	<b>S55</b>	<b>S63</b>	<b>S71</b>	<b>S79</b>	<b>S87</b>
<b>C</b>	<b>PC</b>	<b>S8</b>	<b>S16</b>	<b>S24</b>	<b>S32</b>	<b>S40</b>	<b>S48</b>	<b>S56</b>	<b>S64</b>	<b>S72</b>	<b>S80</b>	<b>S88</b>
<b>D</b>	<b>S1</b>	<b>S9</b>	<b>S17</b>	<b>S25</b>	<b>S33</b>	<b>S41</b>	<b>S49</b>	<b>S57</b>	<b>S65</b>	<b>S73</b>	<b>S81</b>	<b>S89</b>
<b>E</b>	<b>S2</b>	<b>S10</b>	<b>S18</b>	<b>S26</b>	<b>S34</b>	<b>S42</b>	<b>S50</b>	<b>S58</b>	<b>S66</b>	<b>S74</b>	<b>S82</b>	<b>S90</b>
<b>F</b>	<b>S3</b>	<b>S11</b>	<b>S19</b>	<b>S27</b>	<b>S35</b>	<b>S43</b>	<b>S51</b>	<b>S59</b>	<b>S67</b>	<b>S75</b>	<b>S83</b>	<b>S91</b>
<b>G</b>	<b>S4</b>	<b>S12</b>	<b>S20</b>	<b>S28</b>	<b>S36</b>	<b>S44</b>	<b>S52</b>	<b>S60</b>	<b>S68</b>	<b>S76</b>	<b>S84</b>	<b>S92</b>
<b>H</b>	<b>S5</b>	<b>S13</b>	<b>S21</b>	<b>S29</b>	<b>S37</b>	<b>S45</b>	<b>S53</b>	<b>S61</b>	<b>S69</b>	<b>S77</b>	<b>S85</b>	<b>S93</b>

**Buffer** = Sample Diluent

**NC** = Negative Control

**PC** = Positive Control

**S** = Test Samples in singlicate

**C0401D16GB April 2017**

### References

- 1.0 Langedijk J, Schaaper W, Meloen R, van Oirschot J. Proposed three dimensional model for the attachment protein G of respiratory syncytial virus. *Journal of General Virology* 1996, 77, 1249-1257
- 2.0 Quinting B, Beatrice R, Letellier C, Boxus M, Kerkhofs P, Schynts F, Collard A. Development of a 1 step enzyme linked immunosorbent assay for the rapid diagnosis of bovine respiratory syncytial virus in post-mortem specimens. *J Vet Diagn. Invest.* 19: 238-243, 2007.
- 3.0 Verhoeff J, van der Ban M, van Wieuwstadt A. Bovine Respiratory syncytial virus infections in young dairy cattle: clinical and haematological findings. *Vet Rec.* 114: 9-12 1984
- 4.0 Saif, L.J., Bovine Respiratory Coronavirus. *Vet Clin. Health Am Food Anim. Pract.* M2010 July; 26(2) 349-364
- 5.0 Cho K.O. Hoet AE, Loerch SC., Wittiem TE, Saif LJ. Evaluation of concurrent shedding of bovine coronaviruses via respiratory tract and enteric route in feed lot cattle. *Am J Vet Res.* 2001; 62: 1436-41.
- 6.0 Benfield DA, Saif LJ. Cell culture propagation of a coronavirus isolated from cows with winter dysentery. *J. Clin Microbiol.* 1990, 28: 1454-7
- 7.0 Takiuchi E., Fernandes-Barry A., Fernandes-Alfieri A., Filippsen P., Alcindo-Alfieri A. An outbreak of winter dysentery caused by bovine coronavirus in a high production dairy cattle herd from a tropical country. *Braz. Arch. Biol. Technol.* V. 52 n. special: pp. 57-61, Nov. 2009.

## Appendix 1

Spot Number	Antigen
Spot 1	Blank
Spot 2	BCV Haem
Spot 3	BCV N
Spot 4	BCV S1
Spot 5	BRSV-1
Spot 6	BRSV-B
Spot 7	BRSV tF1
Spot 8	BRSV Nuc